

Product sheet

HS-695T | 300211

Description

HS-695T
HS-695T -A, , HMB-45 .

Organism

Tissue

Disease

Metastatic site

Synonyms Hs 695.T, Hs-695-T, Hs 695T, HS 695T, Hs695T, HS695T, Hs695

Age 26

Gender

Ethnicity

Morphology

Growth properties

Citation HS-695T(300211)

Biosafety level 1

NCBI_TaxID 9606

CellosaurusAccession CVCL_0851

Product sheet

HS-695T | 300211

Protein expression	P53				
Isoenzymes	G6PD, B, PGM1, 1, PGM3, 1, ES-D, 1, Me-2, 0, AK-1, 1, GLO-1, 1, : 0.0427				
Tumorigenic	,				
Mutational profile	BRAF V600				
Karyotype	(P19-40) = 52, Y				
Culture Medium	DMEM, w: 4.5g/L , w: 4mM L- , w: 3.7g/L NaHCO3, w: 1.0mM (820300a)				
Supplements	10% FBS				
Dissociation Reagent					
Subculturing	PBS . T25 3~5ml, T75 5~10ml PBS . T25				
Seeding density	2 x 10 ⁴ /cm ²				
Fluid renewal	2~3				
Post-Thaw Recovery	4 5 x 10 ² , 24 .				
Freeze medium	(FBS) + 10% DMSO , CM-1(Cytion 800100)				

Product sheet

HS-695T | 300211

Thawing and
Culturing Cells

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

Incubation
Atmosphere

37°C, 5%_{CO2} .

Flask Coating

Freezing
Procedure

-78°C .

Shipping
Conditions

-78°C .

Storage
Conditions

-150°C -196°C .80°C .

Sterility

PCR .