

MEG-01 | 300482

Description	MEG-01	(CML)	55	1983
	MEG-01			CD41, CD61, CDw14
	MEG-01	,	(PDGF)	GPIIb/IIIa

Organism

Tissue

Disease

Synonyms	Meg-01, MEG01, Meg01
-----------------	----------------------

Age	55
------------	----

Gender

Ethnicity

Morphology

Cell type

Growth properties /

Citation	MEG-01(300482)
-----------------	-----------------

Biosafety level 1

NCBI_TaxID	9606
-------------------	------

CellosaurusAccession CVCL_0425

MEG-01 | 300482

Antigen expression CD41 +, CD61 +, CDw14 +

Culture Medium RPMI 1640, w: 2.0mM , w: 2.0g/L NaHCO3(820700a)

Supplements 10% FBS

Dissociation Reagent

Subculturing 15ml PBS (T25 3~5ml, T75 5~10ml).

Freeze medium (FBS) + 10% DMSO , CM-1(Cytion 800100)

Thawing and Culturing Cells

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

Incubation Atmosphere 37°C, 5% CO₂ .

MEG-01 | 300482

Flask Coating

Freezing Procedure

-78°C

Shipping Conditions

-78°C

Storage Conditions

-150°C -196°C .80°C

/ / HLA

Sterility

PCR