

BNL CL.2 | 305177

Description	BALB/c CDK BNL CL.2	BNL CL.2	,	.
--------------------	---------------------------	----------	---	---

Organism	
Tissue	

Synonyms	BNL-CL.2, BNL CL2, BNL.CL2, BN-CL2, BNCL-2, BNCL2
-----------------	---

Breed/Subspecies	BALB/c
-------------------------	--------

Age	Embryo
------------	--------

Morphology	
-------------------	--

Growth properties	
--------------------------	--

Citation	BNL CL.2(305177)
-----------------	-------------------

Biosafety level	1
------------------------	---

NCBI_TaxID	10090
-------------------	-------

CellosaurusAccession	CVCL_4383
-----------------------------	-----------

Tumorigenic	,	.
--------------------	---	---

BNL CL.2 | 305177

Culture Medium DMEM, w: 4.5g/L , w: 4mM L- , w: 3.7g/L NaHCO₃, w: 1.0mM (820300a)

Supplements 10% FBS

Dissociation Reagent

Subculturing PBS . T25 3~5ml, T75 5~10ml PBS . T25

Fluid renewal 2~3

Freeze medium (FBS) + 10% DMSO , CM-1(Cytion 800100)

- Thawing and Culturing Cells**
1. .
 2. -150°C , 3 .
 3. 37°C 40~60 .
 4. , 70% .
 5. 8ml 15ml .
 6. 300 x g 3 .
 7. 10ml . T25 , T25
 8. .

Incubation Atmosphere 37°C, 5% CO₂ .

Flask Coating

Freezing Procedure -78°C .

BNL CL.2 | 305177

**Shipping
Conditions**

-78°C

**Storage
Conditions**

-150°C -196°C

.80°C

/ / HLA

Sterility

PCR