

Product sheet

P3X63Ag8 | 305171

Description	P3K27	(BALB/c	MOPC-21	)	.	8-	(20	/mL)	HAT
-------------	-------	---------	---------	---	---	----	-----	------	-----

Organism

Tissue	B	,	,
--------	---	---	---

Synonyms P3x63Ag8, P3-x63-Ag8, P3/x63-Ag8, P3/x63 Ag8, P3/x63/Ag8, P3-x63Ag8, P3x63 Ag8, P3x63 Ag8, P3 x 63Ag8, P3 x 63 Ag8, x63-Ag8, x63-AG8, x63-Ag8, P3x63, x63, GM03571

Breed/Subspecies	BALB/c
------------------	--------

Gender

Morphology

Growth properties

Citation	P3x63Ag8 (Cytion	305171)
----------	------------------	---------

Biosafety level 1

NCBI_TaxID	10090
------------	-------

CellosaurusAccession CVCL_3411

Protein expression	,
--------------------	---

Antigen expression H-2d

Product sheet

P3X63Ag8 | 305171

Culture Medium	RPMI 1640, w: 2.0mM , w: 2.0g/L NaHCO3(820700a)
----------------	--

Supplements 10% FBS

Dissociation Reagent	
----------------------	--

Doubling time 16 26

Subculturing	⁵ /ml , , ml .. 1 x 10
--------------	-----------------------------------

Split ratio 2×10⁵ /mL

Seeding density	3~5 × 10 ⁴ /m ²
-----------------	---------------------------------------

Fluid renewal 2~3

Freeze medium	(FBS) + 10% DMSO , CM-1(Cytion 800100)
---------------	---

Thawing and Culturing Cells

1. .
2. -150℃ , 3 .
3. 37℃ 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

Product sheet

P3X63Ag8 | 305171

Incubation Atmosphere	37°C, 5% CO ₂	.
-----------------------	--------------------------	---

Flask Coating

Shipping Conditions	-78°C	.
---------------------	-------	---

Storage Conditions	-150°C	-196°C	.80°C	.
--------------------	--------	--------	-------	---

Sterility	PCR	.
,	.	.