

Product sheet

PBS | 860015

(PBS)

(PBS)	.	pH	.	PBS
PBS	pH	.1	:	
• 8000 mg/L (NaCl) • 200 mg/L (KCl) • 1150 mg/L Na₂HPO₄ (• 200 mg/L KH₂PO₄) (				
	pH	.		
PBS	.	. EDTA	PBS	.
				2
•				
• +2°C ~ +25°C .				
• 2°C~25°C 24 .				
•				
• +2°C ~ +8°C . +37°C .				
• 37°C .				
•				

(mg/L)
200
200
8000
1150