

**B-LCL-HROC59 | 302073**

<b>Description</b>	B-LCL-HROC59	HROC59	B	(TiBc)	-	(EBV)	B	.	
	B-LCL-HROC59	G(IgG)	,	.		B-LCL-HROC59	IgG	TiBc	IgG

**Organism**

**Tissue**

**Disease**

<b>Synonyms</b>	Bc HROC59, TiBcHROC59
-----------------	-----------------------

<b>Age</b>	76
------------	----

**Gender**

**Ethnicity**

**Morphology**

<b>Cell type</b>	B
------------------	---

**Growth properties**

<b>Citation</b>	B-LCL-HROC59(	302073)
-----------------	---------------	---------

<b>Biosafety level</b>	2
------------------------	---

<b>NCBI_TaxID</b>	9606
-------------------	------

<b>CellosaurusAccession</b>	CVCL_A7US
-----------------------------	-----------

**B-LCL-HROC59 | 302073**

**Surface antigens** CD19

**Viruses** : EBV

**Culture Medium** RPMI 1640, w: 2.0mM , w: 2.0g/L NaHCO<sub>3</sub>( 820700a)

**Supplements** 10% FBS

**Subculturing** <sup>5</sup> /ml , , ml .. 1 x 10

**Freeze medium** (FBS ) + 10% DMSO , CM-1(Cytion 800100)

**Thawing and Culturing Cells**

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

**Incubation Atmosphere** 37°C, 5% CO<sub>2</sub> .

**Flask Coating**

**B-LCL-HROC59 | 302073**

**Freezing Procedure**

-78°C

**Shipping Conditions**

-78°C

**Storage Conditions**

-150°C -196°C .80°C

/ / HLA

**Sterility**

PCR

**HLA**

- A\***: '03:01:01, '24:02:01
- B\***: '01:02:01, '27:05:02
- C\***: '02:02:02, '07:02:01
- DRB1\***: '04:01:01, '15:01:01
- DQA1\***: '01:02:01, '03:03:01
- DQB1\***: '03:02:01, '06:02:01
- DPB1\***: '04:01:01, '14:01:01
- E**: '01:03:02