

B-LCL-HROC59 | 302073

Description	B-LCL-HROC59	HROC59	B	(TiBc)	-	(EBV)	B	.	
	B-LCL-HROC59	G(IgG)	,	.		B-LCL-HROC59	IgG	TiBc	IgG

Organism

Tissue

Disease

Synonyms	Bc HROC59, TiBcHROC59
-----------------	-----------------------

Age	76
------------	----

Gender

Ethnicity

Morphology

Cell type	B
------------------	---

Growth properties

Citation	B-LCL-HROC59(302073)
-----------------	---------------	---------

Biosafety level	2
------------------------	---

NCBI_TaxID	9606
-------------------	------

CellosaurusAccession	CVCL_A7US
-----------------------------	-----------

B-LCL-HROC59 | 302073

Surface antigens CD19

Viruses : EBV

Culture Medium RPMI 1640, w: 2.0mM , w: 2.0g/L NaHCO₃(820700a)

Supplements 10% FBS

Subculturing ⁵ /ml , , ml .. 1 x 10

Freeze medium (FBS) + 10% DMSO , CM-1(Cytion 800100)

Thawing and Culturing Cells

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

Incubation Atmosphere 37°C, 5% CO₂' .

B-LCL-HROC59 | 302073

Flask Coating

Freezing Procedure

-78°C

Shipping Conditions

-78°C

Storage Conditions

-150°C -196°C .80°C

/ / HLA

Sterility

PCR

HLA

A*: '03:01:01, '24:02:01
B*: '01:02:01, '27:05:02
C*: '02:02:02, '07:02:01
DRB1*: '04:01:01, '15:01:01
DQA1*: '01:02:01, '03:03:01
DQB1*: '03:02:01, '06:02:01
DPB1*: '04:01:01, '14:01:01
E: '01:03:02