

Product sheet

KYSE-150 | 305087

Description

KYSE-150 (ESCC) .
KYSE-150 , p16(INK4a) . p16 ,
KYSE-150 ESCC .

Organism

Tissue

Disease

Synonyms

KYSE 150, KYSE150, Kyse150, KY150

Age

49

Gender

Ethnicity

Morphology

Growth
properties

Citation

KYSE-150(305087)

Biosafety level

1

NCBI_TaxID

9606

CellosaurusAccession

CVCL_1348

Product sheet

KYSE-150 | 305087

Culture Medium	Ham F12 RPMI 1640 50:50 (820600a 820702a)
----------------	--

Supplements 5% FBS

Dissociation Reagent	
----------------------	--

Doubling time 25

Subculturing	PBS . T25 3~5ml, T75 5~10ml PBS . T25
--------------	---------------------------------------

Fluid renewal 2~3

Freeze medium	(FBS) + 10% DMSO , CM-1(Cytion 800100)
---------------	---

Thawing and Culturing Cells

1. .
2. -150℃ , 3 .
3. 37℃ 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

Incubation Atmosphere	37℃, 5% CO2' .
-----------------------	----------------

Flask Coating

KYSE-150 | 305087

Freezing Procedure	-78°C	.	.
Shipping Conditions	-78°C	.	.
Storage Conditions	-150°C	-196°C	.80°C
Sterility	PCR	.	.