

M2-10B4 | 400428

Products , IV(I(-), VIII(-).

Culture Medium RPMI 1640, w: 2.0mM , w: 2.0g/L NaHCO3(820700a)

Supplements 10% FBS

Dissociation Reagent

Subculturing PBS . T25 3~5ml, T75 5~10ml PBS . T25

Seeding density $1 \times 10^4 / \text{cm}^2$

Fluid renewal 2~3

Post-Thaw Recovery

Freeze medium (FBS) + 10% DMSO , CM-1(Cytion 800100)

M2-10B4 | 400428

Thawing and Culturing Cells	
1.	.
2.	-150°C , 3
3.	37°C 40~60
4.	, 70%
5.	8ml 15ml
6.	300 x g 3
7.	10ml T25 , T25
8.	.

Incubation Atmosphere 37°C, 5% CO₂

Flask Coating	
	.

Freezing Procedure -78°C

Shipping Conditions	
	-78°C

Storage Conditions -150°C -196°C .80°C

/ / HLA

Sterility	
	PCR
	.