

CTLL-2 | 400482

Description	T	-2, CTLL-2	T	.	F4-5	(FLV)	C57BL/6
	CTLL-2	T		.		IL-2	

Organism

Tissue

Synonyms	CTLL 2, CTLL2, CTLL(2)
-----------------	------------------------

Morphology

Cell type

Growth properties

Citation	CTLL-2(400482)
-----------------	---------	---------

Biosafety level	1
------------------------	---

NCBI_TaxID	10090
-------------------	-------

CellosaurusAccession	CVCL_0227
-----------------------------	-----------

Receptors expressed	IL-2
----------------------------	------

Viruses	()	.
----------------	-----	---

Karyotype

CTLL-2 | 400482

Culture Medium i2Cult()

Subculturing , 50% . 5% . . 100 /ml 48

Seeding density 5×10^5 /mL

Fluid renewal 2~3

Post-Thaw Recovery 48 .

Freeze medium (FBS) + 10% DMSO , CM-1(Cytion 800100)

Thawing and Culturing Cells

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

Incubation Atmosphere 37°C, 5% CO₂ .

Flask Coating .

CTLL-2 | 400482

**Freezing
Procedure**

-78°C

**Shipping
Conditions**

-78°C

**Storage
Conditions**

-150°C -196°C .80°C

/ / HLA

Sterility

PCR