

Product sheet

Lec1 | 305010

Description

Lec1 CHO Pro-5 , (wheat germ agglutinin) .
(high-mannose type oligosaccharides) .
Lec1 , N- . Lec1 , ,
(in vitro) .

Organism

Tissue

Synonyms CHO-Lec1, CHO Lec1, Pro-Lec1.3C, Pro-5 Lec1.3c, Pro-5WgaRI3C

Age

Morphology

Growth properties

Citation Lec1 (Cytion 305010)

Biosafety level 1

NCBI_TaxID 10029

CellosaurusAccession CVCL_3440

Culture Medium MEM, w: 2.0mM , : , : , w: 1.0mM , w: 2.2g/L NaHCO3

Product sheet

Lec1 | 305010

Supplements	10% FBS									
Dissociation Reagent										
Subculturing		PBS	. T25	3~5ml, T75	5~10ml	PBS	.	T25		
Seeding density	$2\sim4\times10^4$ /cm ²									
Fluid renewal	2~3									
Freeze medium		(FBS) + 10% DMSO			,		CM-1(Cytion	800100)		
Thawing and Culturing Cells	1.									
	2.	-150°C			,	3				
	3.		37°C			40~60				
	4.			,	70%					
	5.		8ml		15ml					
	6.	300 x g	3							
	7.	10ml		.		T25	,		T25	
	8.									
Incubation Atmosphere	37°C, 5% _{CO2} ' .									
Shipping Conditions		-78°C			.					.
Storage Conditions		-150°C	-196°C		. 80°C					.

Product sheet

Lec1 | 305010

Sterility

PCR