

Product sheet

TRI-C3H | 400359

Description

C3H

Organism

Tissue

Disease

Breed/Subspecies

C3H

Gender

Morphology

Growth
properties

Citation

TRI-C3H(400359)

Biosafety level

1

NCBI_TaxID

10090

CellosaurusAccession

CVCL_5902

Viruses

MAP

Culture
Medium

F12, w: 1.0mM , w: 1.0mM , w: 1.1g/L NaHCO3(820600a)

TRI-C3H | 400359

Supplements 10% FBS

**Dissociation
Reagent**

Subculturing PBS . T25 3~5ml, T75 5~10ml PBS . T25

**Seeding
density** $1 \times 10^4 / \text{cm}^2$

Fluid renewal 1~2

**Post-Thaw
Recovery** $5 \times 10^4 / \text{cm}^2$, 24 .

**Freeze
medium** (FBS) + 10% DMSO , CM-1(Cytion 800100)

**Thawing and
Culturing Cells**

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

**Incubation
Atmosphere** 37°C, 5%_{CO2} .

Flask Coating

Product sheet

TRI-C3H | 400359

Freezing Procedure	-78°C	.	.
--------------------	-------	---	---

Shipping Conditions	-78°C	.	.
---------------------	-------	---	---

Storage Conditions	-150°C	-196°C	.80°C	.
--------------------	--------	--------	-------	---

Sterility	PCR	.
,	.	