

Product sheet

HROC69 T0 M2 | 300829

Description	2006	PD	CRC	.
Organism				
Tissue	, UICC IIIa,	CRC	(	, TNM T3N0M1R0L0V0, G3, Lk(n) + 0, Σ Lk(n) 34)
Disease				
Synonyms	HROC69x			
Age	62			
Gender				
Ethnicity				
Morphology				
Growth properties				
Citation	HROC69 T0 M2(300829)			
Biosafety level	1			
NCBI_TaxID	9606			
CellosaurusAccession	CVCL_1G07			
Protein expression	- , -, (TLR) 3 -, TLR4 , TLR7 , TLR8 -, PTEN			

Product sheet

HROC69 T0 M2 | 300829

Antigen expression	(PS)
Tumorigenic	,
Viruses	SV40, JC/BK, HBV, HCV, HIV .
Ploidy status	
MSI-status	MSS
Mutational profile	APCR1450*, p53R306*, K-Raswt, B-RAFwt, N-Raswt, H-Raswt, SNP rs12628 27, PIK3CAwt
Culture Medium	DMEM: F12(1:1), w: 3.1 g/L , w: 2.5mM L- , w: 15mM HEPES, w: 0.5mM , w: 1.2 g/L NaHCO3 (Cytion
Supplements	10% FBS
Dissociation Reagent	
Doubling time	25~40
Subculturing	PBS . T25 3~5ml, T75 5~10ml PBS . T25
Seeding density	2 x 10 ⁴ /cm ²
Fluid renewal	3~5
Post-Thaw Recovery	2
Freeze medium	(FBS) + 10% DMSO , CM-1(Cytion 800100)

Product sheet

HROC69 T0 M2 | 300829

Thawing and
Culturing Cells

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

Incubation
Atmosphere

37°C, 5%_{CO2} .

Flask Coating

Freezing
Procedure

-78°C .

Shipping
Conditions

-78°C .

Storage
Conditions

-150°C -196°C . 80°C .

Sterility

PCR .