

Product sheet

Colo-320DM | 300153

Description

COLO-320DM 55
COLO-320DM

Organism

Tissue

, C

Disease

Synonyms

COLO_320DM, COLO-320-DM, COLO #320DM, COLO320/DM, COLO320-DM, COLO320DM, COLO320DM, COLO320 DM, COLO 320 DM, COLO 320 (DM), Colorado 320

Age

55

Gender

Ethnicity

Morphology

Growth properties

Citation

COLO-320DM(300153)

Biosafety level

1

NCBI_TaxID

9606

CellosaurusAccession

CVCL_0219

Product sheet

Colo-320DM | 300153

Isoenzymes	PGM1,1, PGM3, 2, G6PD, B, PEP-D, 1, 6PGD, A, ES-D, 1				
Tumorigenic	,				
Products	,	,	,	(ACTH),	
Culture Medium	F12, w: 1.0mM	, w: 1.0mM	, w: 1.1g/L NaHCO3(	820600a)	
Supplements	10% FBS				
Dissociation Reagent					
Subculturing		PBS	. T25	3~5ml, T75	5~10ml PBS . T25
Seeding density	$1 \times 10^4 / \text{cm}^2$				
Fluid renewal	3~5				
Post-Thaw Recovery	$5 \times 10^4 / \text{cm}^2$	,	24	.	
Freeze medium	(FBS) + 10% DMSO , CM-1(Cytion 800100)				

Product sheet

Colo-320DM | 300153

Thawing and
Culturing Cells

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

Incubation
Atmosphere

37°C, 5%_{CO2} .

Flask Coating

Freezing
Procedure

-78°C .

Shipping
Conditions

-78°C .

Storage
Conditions

-150°C -196°C .80°C .

Sterility

PCR .