

(hGF) | 300703

Description

(hGF)
hGF , * *

Organism

Tissue

Gingiva

Applications

,

Cell type

Growth
properties

Citation

(hGF)(300703)

NCBI_TaxID

9606

Culture
Medium

DMEM: F12(1:1), w: 3.1 g/L , w: 2.5mM L- , w: 15mM HEPES, w: 0.5mM , w: 1.2 g/L NaHCO3 (Cytion

Supplements

10% FBS, 10ng/mL bFGF, 10 /L

Dissociation
Reagent

Subculturing

PBS . T25 3~5ml, T75 5~10ml PBS . T25

(hGF) | 300703

Freeze medium	90% FBS + 10% DMSO	,	CM-1(Cytion	800100)
Thawing and Culturing Cells	1.	.		
	2.	-150°C	,	3.
	3.	37°C	40~60	.
	4.	,	70%	.
	5.	8ml	15ml	.
	6.	300 x g	3	.
	7.	10ml	.	T25
	8.		.	T25
Incubation Atmosphere	37°C, 5% CO ₂	.		
Flask Coating				
Freezing Procedure	-78°C	.		.
Shipping Conditions	-78°C	.		.
Storage Conditions	-150°C	-196°C	.80°C	.

(hGF) | 300703

Sterility	PCR
------------------	-----