

HCC78 | 302156

**Description** HCC78 , . HCC78 ROS1 SLC34A2  
HCC78 ROS1 , ROS1

**Organism**

**Tissue**

**Disease**

**Synonyms** HCC-78, HCC0078, 78

**Age** 65

**Gender**

**Ethnicity**

**Growth properties** ,

**Citation** HCC78( 302156)

**Biosafety level** 1

**NCBI\_TaxID** 9606

**CellosaurusAccession** CVCL\_2061

HCC78 | 302156

**Culture Medium** RPMI 1640, w: 2.0mM , w: 2.0g/L NaHCO3( 820700a)

**Supplements** 10% FBS

**Dissociation Reagent**

**Subculturing** PBS . T25 3~5ml, T75 5~10ml PBS . T25

**Freeze medium** (FBS ) + 10% DMSO , CM-1(Cytion 800100)

- Thawing and Culturing Cells**
1. .
  2. -150°C , 3 .
  3. 37°C 40~60 .
  4. , 70% .
  5. 8ml 15ml .
  6. 300 x g 3 .
  7. 10ml . T25 , T25
  8. .

**Incubation Atmosphere** 37°C, 5% CO2' .

**Flask Coating** .

**Freezing Procedure** -78°C .

HCC78 | 302156

---

**Shipping  
Conditions**

-78°C

**Storage  
Conditions**

-150°C -196°C .80°C

/ / HLA

**Sterility**

PCR