

TK6 | 300357

Description

TK6 5 (TK)
tk 3 TK6

Organism

Tissue

Synonyms TK-6, H2BT

Age 5

Gender

Cell type

Growth properties

Citation TK6(300357)

Biosafety level 2

NCBI_TaxID 9606

CellosaurusAccession CVCL_0561

Culture Medium RPMI 1640, w: 2.0mM , w: 2.0g/L NaHCO3(820700a)

TK6 | 300357

Supplements	10 % FBS, 2,5 %		
Subculturing	5×10^5 /mL	5×10^6	$1 \times 10^5 \sim 1 \times 10^6$ /mL
Seeding density	1×10^5 /mL		
Fluid renewal	2~3		
Freeze medium	(FBS) + 10% DMSO		CM-1(Cytion 800100)
Thawing and Culturing Cells	<ol style="list-style-type: none"> 1. 2. -150°C 3. 37°C 40~60 4. 70% 5. 8ml 15ml 6. 300 x g 3 7. 10ml T25 T25 8. 		
Incubation Atmosphere	37°C, 5% CO ₂		
Flask Coating			
Freezing Procedure	-78°C		
Shipping Conditions	-78°C		

