



**RBL-2H3 | 305194**

**Culture Medium** EMEM(MEM ), w: 2mM L- , w: 2.2g/L NaHCO3, w: EBSS( 820100a)

**Supplements** 10% FBS 1% NEAA

**Dissociation Reagent**

**Subculturing** PBS . T25 3~5ml, T75 5~10ml PBS . T25

**Split ratio** 1:2 ~ 1:4

**Fluid renewal** 2~3

**Freeze medium** (FBS ) + 10% DMSO , CM-1(Cytion 800100)

- Thawing and Culturing Cells**
1. .
  2. -150°C , 3 .
  3. 37°C 40~60 .
  4. , 70% .
  5. 8ml 15ml .
  6. 300 x g 3 .
  7. 10ml . T25 , T25
  8. .

**Incubation Atmosphere** 37°C, 5% CO2' .

**Flask Coating**

RBL-2H3 | 305194

<b>Freezing Procedure</b>	-78°C	.
---------------------------	-------	---

<b>Shipping Conditions</b>	-78°C	.
----------------------------	-------	---

<b>Storage Conditions</b>	-150°C	-196°C	.80°C	.
---------------------------	--------	--------	-------	---

/ / HLA

<b>Sterility</b>	PCR	.
------------------	-----	---