

CEM/C1 | 305103

Description	CEM/C1	,	II	CCRF-CEM	T	.
	CEM/C1	T	T	.	,	,
	CEM/C1			.		

Organism

Tissue	
---------------	--

Disease	T
----------------	---

Synonyms	CCRF-CEM C1, CEM-C1, CEM.C1, CEMC1
-----------------	------------------------------------

Age	4
------------	---

Gender

Morphology	
-------------------	--

Growth properties

Citation	CEM/C1(305103)
-----------------	---------	---------

Biosafety level	1
------------------------	---

NCBI_TaxID	9606
-------------------	------

CellosaurusAccession	CVCL_3496
-----------------------------	-----------

CEM/C1 | 305103

Culture Medium	RPMI 1640, w: 2.0mM	, w: 2.0g/L NaHCO ₃ (820700a)
Supplements	10%	FBS	
Subculturing	⁵ /ml	, ,	ml .. 1 x 10
Fluid renewal	2~3		
Freeze medium		(FBS) + 10% DMSO	, CM-1(Cytion 800100)
Thawing and Culturing Cells	1.		.
	2.	-150°C	, 3 .
	3.	37°C	40~60 .
	4.	, 70%	.
	5.	8ml	15ml .
	6.	300 x g 3	.
	7.	10ml	. T25 , T25
	8.		.
Incubation Atmosphere	37°C, 5% CO ₂	.	
Flask Coating			
Freezing Procedure		-78°C	.
Shipping Conditions		-78°C	.

Product sheet

CEM/C1 | 305103

Storage Conditions	-150°C	-196°C	-80°C
---------------------------	--------	--------	-------

/ / HLA

Sterility	PCR
------------------	-----