

iPSC-hDPSC | 300622

Description

IPSC-hDPSC (hDPSC) . 70~80%

HDPSC (iPSC) . (Oct-4, Sox-2, Klf-4) p53 , EBNA

Organism

Tissue

Disease

Applications , , , ,

Ethnicity

Morphology IPSC,

Cell type

Growth properties

Citation iPSC-hDPSC(300622)

Biosafety level 1

NCBI_TaxID 9606

CellosaurusAccession CVCL_B7DW

GMO Status GMO-S1: iPSC (iPSC-hDPSC) OCT4, SOX2, KLF4 c-MYC

3 | 69124 | |
: +49(0)6221 405780 | www.cytion.com | info@cytion.com

Product sheet

iPSC-hDPSC | 300622

Thawing and
Culturing Cells

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

Incubation
Atmosphere

37°C, 5% CO₂ , .

Flask Coating

Freezing
Procedure

-78°C . .

Shipping
Conditions

-78°C . .

Storage
Conditions

-150°C -196°C . 80°C .

Sterility

PCR .
 , .