



PIEC | 305213

<b>Subculturing</b>		PBS	T25	3~5ml, T75	5~10ml	PBS	T25
<b>Split ratio</b>	1:2 ~ 1:4						
<b>Fluid renewal</b>	2~3						
<b>Freeze medium</b>		(FBS ) + 10% DMSO			CM-1(		800100)
<b>Thawing and Culturing Cells</b>	1.						
	2.	-150°C			3		
	3.		37°C		40~60		
	4.			70%			
	5.		8ml		15ml		
	6.	300 x g	3				
	7.		10ml		T25		T25
	8.						
<b>Incubation Atmosphere</b>	37°C, 5% CO <sub>2</sub>						
<b>Flask Coating</b>							
<b>Freezing Procedure</b>		-78°C					
<b>Shipping Conditions</b>		-78°C					

PIEC | 305213

**Storage  
Conditions**

-150°C -196°C .80°C

/ / HLA

**Sterility**

PCR