

HROC278Met1 T2 M2 | 300836

Description 2006 PD CRC .

Organism

Tissue , CRC (, TNM T4N2M1R0L1V1, G3, Lk(n) +19, Σ Lk(n) 29)

Disease

Age 76

Gender

Ethnicity

Morphology

Growth properties

Citation HROC278Met1 T2 M2(300836)

Biosafety level 1

NCBI_TaxID 9606

CellosaurusAccession CVCL_1U90

Protein expression PTEN

Tumorigenic ,

HROC278Met1 T2 M2 | 300836

Viruses SV40, JC/BK, HBV, HCV, HIV .

MSI-status MSI-L

Mutational profile B-RAFV600E APCwt, p53wt, K- , N- , H- , PIK3CAwt

Culture Medium DMEM: F12(1:1), w: 3.1 g/L , w: 2.5mM L- , w: 15mM HEPES, w: 0.5mM , w: 1.2 g/L NaHCO3 (Cytion

Supplements 10% FBS

Dissociation Reagent

Doubling time 29

Subculturing PBS . T25 3~5ml, T75 5~10ml PBS . T25

Seeding density 2×10^4 /cm²

Fluid renewal 3~5

Post-Thaw Recovery

Freeze medium (FBS) + 10% DMSO , CM-1(Cytion 800100)

HROC278Met1 T2 M2 | 300836

Thawing and Culturing Cells

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

Incubation Atmosphere

37°C, 5% CO₂ .

Flask Coating

.

Freezing Procedure

-78°C .

Shipping Conditions

-78°C .

Storage Conditions

-150°C -196°C .80°C .

/ / HLA

Sterility

PCR .

,

.