

NCI-H2452 | 300391

**Description**

NCI-H2452

NCI-H2452

(CAR) CD46

NCI-H2452

NCI-H2452

(PGRN)

**Organism**

**Tissue**

**Disease**

**Synonyms**

NCI-H2452, H-2452, NCIH2452

**Age**

**Gender**

**Ethnicity**

**Morphology**

**Growth properties**

**Citation**

NCI-H2452( 300391)

**NCBI\_TaxID**

9606

**CellosaurusAccession**

CVCL\_1553

NCI-H2452 | 300391

**Culture Medium** RPMI 1640, w: 2.0mM , w: 2.0g/L NaHCO3( 820700a)

**Supplements** 10% FBS

**Dissociation Reagent**

**Subculturing** PBS . T25 3~5ml, T75 5~10ml PBS . T25

**Freeze medium** (FBS ) + 10% DMSO , CM-1(Cytion 800100)

- Thawing and Culturing Cells**
1. .
  2. -150°C , 3 .
  3. 37°C 40~60 .
  4. , 70% .
  5. 8ml 15ml .
  6. 300 x g 3 .
  7. 10ml . T25 , T25
  8. .

**Incubation Atmosphere** 37°C, 5% CO2' .

**Flask Coating** .

**Freezing Procedure** -78°C .

NCI-H2452 | 300391

---

**Shipping  
Conditions**

-78°C

**Storage  
Conditions**

-150°C

-196°C

-80°C

/

/ HLA

**Sterility**

PCR