

**PATU8988T | 305133**

**Description** 1985 64 . PATU8988T PATU8988S(DSM ACC 204)

**Organism**

**Tissue**

**Disease**

**Synonyms** PA-TU-8988T, PaTu8988T, PaTu8988T, PATU8988T, PaTu 8988 T, PaTu-8988t, PATU-8988T, PATU-T, PA-TU T, 8988T, PaCL4

**Age** 64

**Gender**

**Ethnicity**

**Morphology**

**Growth properties**

**Citation** PATU8988T( 305133)

**Biosafety level** 1

**NCBI\_TaxID** 9606

**CellosaurusAccession** CVCL\_1847

**PATU8988T | 305133**

**Culture Medium** DMEM, w: 4.5g/L , w: 4mM L- , w: 3.7g/L NaHCO3, w: 1.0mM ( 820300a)

**Supplements** 10% FBS

**Dissociation Reagent**

**Doubling time** 22~30

**Subculturing** PBS . T25 3~5ml, T75 5~10ml PBS . T25

**Fluid renewal** 2~3

**Freeze medium** (FBS ) + 10% DMSO , CM-1(Cytion 800100)

- Thawing and Culturing Cells**
1. .
  2. -150°C , 3 .
  3. 37°C 40~60 .
  4. , 70% .
  5. 8ml 15ml .
  6. 300 x g 3 .
  7. 10ml . T25 , T25
  8. .

**Incubation Atmosphere** 37°C, 5% CO2' .

**Flask Coating**

PATU8988T | 305133

|                           |       |   |
|---------------------------|-------|---|
| <b>Freezing Procedure</b> | -78°C | . |
|---------------------------|-------|---|

|                            |       |   |
|----------------------------|-------|---|
| <b>Shipping Conditions</b> | -78°C | . |
|----------------------------|-------|---|

|                           |        |        |       |   |
|---------------------------|--------|--------|-------|---|
| <b>Storage Conditions</b> | -150°C | -196°C | .80°C | . |
|---------------------------|--------|--------|-------|---|

/ / HLA

|                  |     |   |
|------------------|-----|---|
| <b>Sterility</b> | PCR | . |
|------------------|-----|---|