

| 305071

|                    |          |   |                  |               |                  |   |
|--------------------|----------|---|------------------|---------------|------------------|---|
| <b>Description</b> | Farage   | B | B                |               |                  |   |
|                    | -4(IL-4) |   | CD23, CD54, CD58 |               | CD21, CD22, CD38 |   |
|                    |          |   | (TdT)            | (RAG-1 RAG-2) | B                | B |

|                 |  |
|-----------------|--|
| <b>Organism</b> |  |
| <b>Tissue</b>   |  |

|                        |   |   |
|------------------------|---|---|
| <b>Disease</b>         | B | B |
| <b>Metastatic site</b> |   |   |

|                 |       |
|-----------------|-------|
| <b>Synonyms</b> | , OL, |
|-----------------|-------|

|            |    |
|------------|----|
| <b>Age</b> | 70 |
|------------|----|

|                  |  |
|------------------|--|
| <b>Gender</b>    |  |
| <b>Ethnicity</b> |  |

|                          |  |
|--------------------------|--|
| <b>Morphology</b>        |  |
| <b>Growth properties</b> |  |

|                 |                 |
|-----------------|-----------------|
| <b>Citation</b> | Farage( 305071) |
|-----------------|-----------------|

|                        |   |
|------------------------|---|
| <b>Biosafety level</b> | 2 |
|------------------------|---|

|                   |      |
|-------------------|------|
| <b>NCBI_TaxID</b> | 9606 |
|-------------------|------|

|                             |           |
|-----------------------------|-----------|
| <b>CellosaurusAccession</b> | CVCL_3302 |
|-----------------------------|-----------|

| 305071

**Culture Medium** RPMI 1640, w: 2.0mM , w: 2.0g/L NaHCO3( 820700a)

**Supplements** 10% FBS 2.5g/L 10mM HEPES

**Doubling time** 48

**Subculturing** 1.5-2 x 10<sup>6</sup> /ml , ml

**Seeding density** 5 x 10<sup>5</sup> /mL

**Fluid renewal** 2~3

**Freeze medium** (FBS ) + 10% DMSO , CM-1(Cytion 800100)

- Thawing and Culturing Cells**
1. .
  2. -150°C , 3 .
  3. 37°C 40~60 .
  4. , 70% .
  5. 8ml 15ml .
  6. 300 x g 3 .
  7. 10ml . T25 , T25
  8. .

**Incubation Atmosphere** 37°C, 5% CO<sub>2</sub>' .

| 305071

**Flask Coating**

**Freezing Procedure**

-78°C

**Shipping Conditions**

-78°C

**Storage Conditions**

-150°C -196°C .80°C

/ / HLA

**Sterility**

PCR