

EB3 | 300373

Description

EB3

EB3

Organism

Tissue

Bone

Disease

Metastatic site

Bone

Applications

3D ,

Synonyms

EB-3, - -3, GM04679

Age

3

Gender

Ethnicity

Morphology

Cell type

B

Growth properties

Citation

EB3(300373)

Biosafety level

2

NCBI_TaxID

9606

EB3 | 300373

CellosaurusAccession CVCL_1185

Surface antigens HLA A3, Aw32, Cw2

Isoenzymes G6PD, A

Viruses EBV(EBNA)

Culture Medium RPMI 1640, w: 2.0mM , w: 2.0g/L NaHCO3(820700a)

Supplements 10% FBS

Subculturing ⁵ /ml , , ml .. 1 x 10

Freeze medium (FBS) + 10% DMSO , CM-1(Cytion 800100)

Thawing and Culturing Cells

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

EB3 | 300373

Incubation Atmosphere 37°C, 5% CO₂

Flask Coating

Freezing Procedure -78°C

Shipping Conditions -78°C

Storage Conditions -150°C -196°C -80°C

/ / HLA

Sterility PCR