

OKT 3 | 300499

Description P3x63Ag8.U1 Muromonab-CD3

Organism

Tissue

Applications

Synonyms OKT3

Breed/Subspecies CAF1 x BALB/c

Morphology

Cell type (, B)

Growth properties /

Citation OKT 3(300499)

Biosafety level 1

NCBI_TaxID 10090

CellSaurusAccession CVCL_2665

Products , T (T), CD3 ,

OKT 3 | 300499

Culture Medium RPMI 1640, w: 2.0mM , w: 2.0g/L NaHCO3(820700a)

Supplements 10% FBS

Subculturing ⁵ /ml , , ml .. 1 x 10

Freeze medium (FBS) + 10% DMSO , CM-1(Cytion 800100)

Thawing and Culturing Cells

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

Incubation Atmosphere 37°C, 5% CO₂ .

Flask Coating .

Freezing Procedure -78°C .

Shipping Conditions -78°C .

Product sheet

OKT 3 | 300499

**Storage
Conditions**

-150°C -196°C .80°C

/ / **HLA**

Sterility

PCR