

PtK1 | 608393

Description	PtK1	PtK1
	PtK1	PtK1

Organism

Tissue

Synonyms Pt K1(NBL-3), NBL-3, PTK-1, PTK 1, PtK 1, PTK 1, PtK1, Pt-K1, Ptk1, 1

Age

Gender

Morphology

Growth properties

Citation PtK1(608393)

Biosafety level 1

NCBI_TaxID 9310

CellosaurusAccession CVCL_0489

Virus susceptibility ()

PtK1 | 608393

Virus resistance 2

Reverse transcriptase

Products

Culture Medium RPMI 1640, w: 2.0mM , w: 2.0g/L NaHCO3(820700a)

Supplements 10% FBS

Dissociation Reagent

Subculturing PBS . T25 3~5ml, T75 5~10ml PBS . T25

Fluid renewal 2

Freeze medium (FBS) + 10% DMSO , CM-1(Cytion 800100)

PtK1 | 608393

Thawing and Culturing Cells	
1.	.
2.	-150°C , 3 .
3.	37°C 40~60 .
4.	, 70% .
5.	8ml 15ml .
6.	300 x g 3 .
7.	10ml . T25 , T25
8.	.

Incubation Atmosphere

37°C, 5% CO₂

Flask Coating	
	.

Freezing Procedure

-78°C

Shipping Conditions	
	-78°C .

Storage Conditions

-150°C -196°C .80°C

/ / HLA

Sterility	
	PCR .
	.