

TM3 | 305167

Description TM3 11~13 , .

Organism

Tissue

Synonyms TM-3

Breed/Subspecies BALB/c

Age 11~13

Gender

Morphology

Growth properties

Citation TM3(305167)

Biosafety level 1

NCBI_TaxID 10090

CellSaurusAccession CVCL_4326

Culture Medium DMEM: F12(1:1), w: 3.1 g/L , w: 2.5mM L- , w: 15mM HEPES, w: 0.5mM , w: 1.2 g/L NaHCO3 (Cytion

TM3 | 305167

Supplements 2.5 % FBS, 5 %

Dissociation Reagent

Subculturing PBS . T25 3~5ml, T75 5~10ml PBS . T25

Fluid renewal 2~3

Freeze medium (FBS) + 10% DMSO , CM-1(Cytion 800100)

Thawing and Culturing Cells

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

Incubation Atmosphere 37°C, 5% CO₂ .

Flask Coating

Freezing Procedure -78°C .

Shipping Conditions -78°C .

TM3 | 305167

**Storage
Conditions**

-150°C -196°C .80°C

/ / HLA

Sterility

PCR