

B95-8 | 601102

<b>Description</b>	B95-8	B	,	(Saguinus oedipus)	.	B
	B95-8		-	.		EBV

**Organism**

**Tissue**

**Synonyms** B95.8, B 95.8, B 95-8, B-95-8, B958, GM07404, GM07404A, GM07404d

**Gender**

**Morphology**

**Growth properties**

**Citation** B95-8( 601102)

**Biosafety level** 2

**NCBI\_TaxID** 9490

**CellosaurusAccession** CVCL\_1953

**Culture Medium** RPMI 1640, w: 2.0mM , w: 2.0g/L NaHCO3( 820700a)

**Supplements** 10% FBS

**B95-8 | 601102**

**Subculturing** 5 /ml , , ml .. 1 x 10

**Fluid renewal** 2~3

**Freeze medium** (FBS ) + 10% DMSO , CM-1(Cytion 800100)

**Thawing and Culturing Cells**

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

**Incubation Atmosphere** 37°C, 5% CO<sub>2</sub> .

**Flask Coating**

**Freezing Procedure** -78°C .

**Shipping Conditions** -78°C .

**Storage Conditions** -150°C -196°C . 80°C .

/ / HLA

B95-8 | 601102

**Sterility**

PCR