

LM/TK(LMTK-) | 305176

Description LM/TK-(LMTK-) (TK) , DNA
HSV-TK HSV-TK

Organism

Tissue ,

Synonyms L-M[TK-], LM TK , L-M(TK-), L M(TK-), LM(TK-), LM(tk-), LM-TK-, LMTK-, L (TK-), L(TK-), L(TK-), L(tk-)

Breed/Subspecies C3H/An

Age 100

Gender

Morphology

Growth properties

Citation LM/TK(LMTK-)(305176)

Biosafety level 1

NCBI_TaxID 10090

CellosaurusAccession CVCL_4536

Antigen expression H-2k

LM/TK(LMTK-) | 305176

Tumorigenic , (1×10⁷ 21 100% (5/5)).

Culture Medium DMEM, w: 4.5g/L , w: 4mM L- , w: 3.7g/L NaHCO₃, w: 1.0mM (820300a)

Supplements 10% FBS

Dissociation Reagent

Subculturing PBS . T25 3~5ml, T75 5~10ml PBS . T25

Fluid renewal 2

Freeze medium (FBS) + 10% DMSO , CM-1(Cytion 800100)

Thawing and Culturing Cells

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

Incubation Atmosphere 37°C, 5% CO₂' .

Flask Coating

LM/TK(LMTK-) | 305176

**Freezing
Procedure**

-78°C

**Shipping
Conditions**

-78°C

**Storage
Conditions**

-150°C -196°C .80°C

/ / HLA

Sterility

PCR