

HROC32 T3 M1 | 300819

Description 2006 PD CRC HROC32

Organism

Tissue , UICC IV, CRC (, TNM T4N2M1R0L0V1 G2, Lk(n) + 9, Σ Lk(n) 14)

Disease

Synonyms HROC32x

Age 82

Gender

Ethnicity

Morphology

Growth properties

Citation HROC32 T3 M1(300819)

Biosafety level 1

NCBI_TaxID 9606

CellosaurusAccession CVCL_1D07

Protein expression PTEN

HROC32 T3 M1 | 300819

Antigen expression	CD15 +, CD24 +, CD44 +, CD55 +, CD58 +, CD50 +, CD 54 +, CD66acde +, CD71 +, CD102 +, CD326 +, CD80 -, CD86-, EpCAM +, HLA-A2 +		
Tumorigenic	,		
Viruses	SV40, JC/BK, HBV, HCV, HIV .		
Ploidy status			
MSI-status	MSS		
Mutational profile	APCwt, p53R282W, K-RasG12A, N-Raswt, H-Raswt 27 SNP rs12628, PIK3CAst, BRafwt		
Culture Medium	DMEM: F12(1:1), w: 3.1 g/L , w: 2.5mM L- , w: 15mM HEPES, w: 0.5mM , w: 1.2 g/L NaHCO3 (Cytion		
Supplements	10% FBS		
Dissociation Reagent			
Doubling time	30		
Subculturing	PBS	. T25	3~5ml, T75 5~10ml PBS . T25
Seeding density	2 x 10 ⁴ /cm ²		
Fluid renewal	3~5		
Post-Thaw Recovery	1~2		
Freeze medium	(FBS) + 10% DMSO , CM-1(Cytion 800100)		

HROC32 T3 M1 | 300819

Thawing and Culturing Cells

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

Incubation Atmosphere 37°C, 5% CO₂ .

Flask Coating

.

Freezing Procedure -78°C .

Shipping Conditions

-78°C .

Storage Conditions -150°C -196°C .80°C .

/ / HLA

Sterility PCR .