

Product sheet

RPMI 1640, w: 2.0mM , w: 2.0g/L NaHCO3 | 820700a

RPMI RPMI 1640 .1966 E. , E. ,H.
RPMI 1640 RPMI 1640 HeLa, Jurkat, MCF-7, PC12, PB
RPMI 1640 , , B12, PABA
RPMI 1640 pH 5-10% CO2 .

•

- +2°C ~ +8°C .
- 4°C 6~8 .

•

- +2°C~+8°C . +37°C .
- 37°C .
- .

(mg/L)		
10.00		
L- -L-	434.40	
L-	200.00	
L- H2O	56.82	
L-	20.00	
L- 2HCl	65.20	
L-	20.00	
L- HCl H2O	20.27	
L- -L-	20.00	

Product sheet

RPMI 1640, w: 2.0mM , w: 2.0g/L NaHCO3 | 820700a

L-	50.00
L-	50.00
L- HCl	40.00
L-	15.00
L-	15.00
L-	20.00
L-	30.00
L-	20.00
L-	5.00
L- $2Na_2H_2O$	28.83
L-	20.00
p-	1.00
D-	0.20
	3.00
	0.25
	1.00
	35.00
	1.00
HCl	1.00
	0.20
HCl	1.00

Product sheet

RPMI 1640, w: 2.0mM , w: 2.0g/L NaHCO3 | 820700a

B12	0.005
Ca(NO3)2 4H2O	100.00
KCl	400.00
MgSO4 7H2O	100.00
NaCl	6000.00
NaHCO3	2000.00
Na2HPO4	800.00
D-	2000.00
L-	1.00
	5.30