

D341Med | 305136

<b>Description</b>	D341 Med	1988 Friedman	3
		D341 Med	

**Organism**

<b>Tissue</b>	,
---------------	---

**Disease**

<b>Synonyms</b>	D-341 , D-341 , D-341 , D341_ , D341 , D341 , D341 , D-341, D341, 341, H341
-----------------	---

<b>Age</b>	3,5
------------	-----

**Gender**

**Ethnicity**

**Morphology**

**Growth properties**

<b>Citation</b>	D341Med( 305136)
-----------------	------------------

<b>Biosafety level</b>	1
------------------------	---

<b>NCBI_TaxID</b>	9606
-------------------	------

<b>CellosaurusAccession</b>	CVCL_0018
-----------------------------	-----------

**D341Med | 305136**

**Protein expression** , , , S100(S-100) , , UJ13A

**Tumorigenic**

**Culture Medium** EMEM(MEM ), w: 2mM L- , w: 2.2g/L NaHCO<sub>3</sub>, w: EBSS( 820100a)

**Supplements** 10% FBS 1% NEAA

**Doubling time** 37

**Subculturing** <sup>5</sup> /ml , , ml .. 1 x 10

**Fluid renewal** 2~3

**Freeze medium** (FBS ) + 10% DMSO , CM-1(Cytion 800100)

**Thawing and Culturing Cells**

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

**Incubation Atmosphere** 37°C, 5% CO<sub>2</sub> .

D341Med | 305136

**Flask Coating**

**Freezing Procedure**

-78°C

**Shipping Conditions**

-78°C

**Storage Conditions**

-150°C -196°C .80°C

/ / HLA

**Sterility**

PCR