

U2OS-CRISPR-SNAPf-Nup133 | 300666

**Description** U2OS-CRISPR-SNAPf-Nup133 U2OS , NUP133 CRISPR/Cas9 C- SN  
SNAPf NPC O6- -DNA SNAP- SNAP  
NPC NPC , Y- . U2OS

**Organism**

**Tissue** Bone

**Disease**

**Age** 15

**Gender**

**Ethnicity**

**Morphology**

**Growth properties**

**Citation** U2OS-CRISPR-SNAPf-Nup133( 300666)

**Biosafety level** 1

**NCBI\_TaxID** 9606

**Depositor** (EMBL)

**GMO Status** GMO-S1: (U2OS-CRISPR-SNAPf-Nup133) CRISPR SNAPf-Nup133 Nup133

**U2OS-CRISPR-SNAPf-Nup133 | 300666**

**Protein expression** Nup133, SNAPf-tag

**Culture Medium** 5a, w: 3.0g/L , w: , w: 2.0mM , w: 2.2g/L NaHCO3( 820200a)

**Supplements** 10% FBS, 3.0g/L , , 2.0mM , 2.2g/L NaHCO3, 1% NEAA

**Dissociation Reagent**

**Subculturing** PBS . T25 3~5ml, T75 5~10ml PBS . T25

**Freeze medium** (FBS ) + 10% DMSO , CM-1(Cytion 800100)

- Thawing and Culturing Cells**
1. .
  2. -150°C , 3 .
  3. 37°C 40~60 .
  4. , 70% .
  5. 8ml 15ml .
  6. 300 x g 3 .
  7. 10ml . T25 , T25
  8. .

**Incubation Atmosphere** 37°C, 5% CO2' .

U2OS-CRISPR-SNAPf-Nup133 | 300666

Flask Coating

Freezing Procedure

-78°C

Shipping Conditions

-78°C

Storage Conditions

-150°C -196°C .80°C

/ / HLA

Sterility

PCR