

- | 300645

Description

(hMSCs) , , . (stromal vasculature)
hMSC , . , ,
hMSC , . , , ,
hMSC hMSC .

Organism

Tissue

Applications

, ,

Age

Gender

Ethnicity

Morphology

5 , . 2% .

Cell type

Growth properties

Citation

, (300645)

Biosafety level 1

NCBI_TaxID

9606

Product sheet

- | 300645

Antigen expression	CD73/CD90/CD105() CD14/CD34/CD45/HLA-DR()		MSC(P2-P3)	
Viruses	HBV(PCR),	(PCR), HIV-1/2(IFA)	HBV, HCV, HSV1, HSV2, CMV, EBV, HHV6,	
Culture Medium	MEM, w: 2.0mM	, :	, w: 1.0mM	, w: 2.2g/L NaHCO3
Supplements	10% FBS, 2ng/mL bFGF			
Dissociation Reagent	-EDTA			
Subculturing	:	PBS	. PBS	CO2 /EDTA 2-3 T25
Seeding density	$1\sim3 \times 10^4 /cm^2$			
Fluid renewal	24	2~3	.	
Freeze medium	80% FBS + 10%	+ 10% DMSO	,	CM-1(800100)

- | 300645

**Thawing and
Culturing Cells**

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

**Incubation
Atmosphere**

37°C, 5%_{CO2} .

Flask Coating

**Freezing
Procedure**

-78°C .

**Shipping
Conditions**

-78°C .

**Storage
Conditions**

-150°C -196°C . 80°C .

Sterility

PCR .