

IEC-6 | 302149

Description IEC-6, . IEC-6
IEC-6 (IBD) .

Organism

Tissue

Applications .

Synonyms IEC 6, IEC6, 6

Breed/Subspecies (CD(SD))

Age 18-24

Gender

Morphology

Cell type

Growth properties

Citation IEC-6(302149)

Biosafety level 1

NCBI_TaxID 10116

CellosaurusAccession CVCL_0343

IEC-6 | 302149

Culture Medium DMEM, w: 4.5g/L , w: 4mM L- , w: 3.7g/L NaHCO₃, w: 1.0mM (820300a)

Supplements 10% FBS

Dissociation Reagent

Subculturing PBS . T25 3~5ml, T75 5~10ml PBS . T25

Freeze medium (FBS) + 10% DMSO , CM-1(Cytion 800100)

Thawing and Culturing Cells

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

Incubation Atmosphere 37°C, 5% CO₂ .

Flask Coating

IEC-6 | 302149

Freezing Procedure	-78°C	.
---------------------------	-------	---

Shipping Conditions	-78°C	.
----------------------------	-------	---

Storage Conditions	-150°C	-196°C	.80°C	.
---------------------------	--------	--------	-------	---

/ / HLA

Sterility	PCR	.
------------------	-----	---