

HROC43 | 300823

Description 2006 PD CRC .

Organism

Tissue , UICC IIIb

Disease , TNM T3N2M0R0L1V0, G3, Lk(n) + 30, Σ Lk(n) 36

Age 72

Gender

Ethnicity

Morphology

Growth properties

Citation HROC43(300823)

Biosafety level 1

NCBI_TaxID 9606

CellosaurusAccession CVCL_1D20

Protein expression PTEN

Antigen expression CD326+, CD44+, CD15+, CD71+, CD73+, CD274+, CD47+, CD54+, CD95+, CD276+, CD133- , CD66acdeweak, IDO+, cFLIP+, MHC-I+, MHCII IFN-γ , EpCAM+

HROC43 | 300823

Tumorigenic ,

Viruses SV40, JC/BK, HBV, HCV, HIV .

Ploidy status

MSI-status MSS

Mutational profile APCQ1429*, p53S241fs*5, K- , N- , H- , PIK3CAwt, B-

Culture Medium DMEM: F12(1:1), w: 3.1 g/L , w: 2.5mM L- , w: 15mM HEPES, w: 0.5mM , w: 1.2 g/L NaHCO3 (Cytion 820

Supplements 10% FBS

Dissociation Reagent

Doubling time 27

Subculturing PBS . T25 3~5ml, T75 5~10ml PBS . T25

Seeding density 2×10^4 /cm²

Fluid renewal 3~5

Post-Thaw Recovery 1

Freeze medium (FBS) + 10% DMSO , CM-1(Cytion 800100)

HROC43 | 300823

Thawing and Culturing Cells

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

Incubation Atmosphere 37°C, 5%_{CO2} .

Flask Coating

Freezing Procedure -78°C .

Shipping Conditions

-78°C .

Storage Conditions -150°C -196°C .80°C .

/ / HLA

Sterility PCR .