

B-LCL-HROC43 | 302067

Description	B-LCL-HROC43	B	-	(EBV)	B	.	A
	, B-LCL-HROC43	CD19	CD20	CD23, CD80		B	. MHC
	, B-LCL-HROC43	(: IgG, IgM	IgA)	,	.		

Organism	
Tissue	

Disease	
Synonyms	Bc HROC43

Age	72
------------	----

Gender	
Ethnicity	

Morphology	
Cell type	B

Growth properties	
--------------------------	--

Citation	B-LCL-HROC43(302067)
-----------------	-----------------------

Biosafety level	2
------------------------	---

NCBI_TaxID	9606
-------------------	------

CellosaurusAccession	CVCL_A7UN
-----------------------------	-----------

B-LCL-HROC43 | 302067

Surface antigens CD19

Viruses : EBV

Culture Medium RPMI 1640, w: 2.0mM , w: 2.0g/L NaHCO₃(820700a)

Supplements 10% FBS

Subculturing ⁵ /ml , , ml .. 1 x 10

Freeze medium (FBS) + 10% DMSO , CM-1(Cytion 800100)

Thawing and Culturing Cells

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

Incubation Atmosphere 37°C, 5% CO₂ .

Flask Coating

B-LCL-HROC43 | 302067

Freezing Procedure

-78°C

Shipping Conditions

-78°C

Storage Conditions

-150°C -196°C .80°C

/ / HLA

Sterility

PCR

HLA

- A***: '01:01:01, '02:01:01
- B***: '15:01:01, '37:01:01
- C***: '03:03:01, '06:02:01
- DRB1***: '10:01:01, '12:01:01
- DQA1***: '01:05:01, '05:05:01
- DQB1***: '03:01:01, '05:01:01
- DPB1***: '01:01:01, '05:01:01
- E**: '01:01:01