

Hep-64.1 | 400205

Description Hep-64.1 , C57BL/6Jmouse . Hep-64.1
Hep-64.1 . Hep-64.1 DNA
p53 . Hep-64.1 p53 c-Ha-ras

Organism

Tissue

Disease

Synonyms HEP-64.1, 64.1

Breed/Subspecies C57BL/6J

Age

Gender

Morphology

Growth properties

Citation Hep-64.1(400205)

Biosafety level 1

NCBI_TaxID 10090

CellosaurusAccession CVCL_5770

Hep-64.1 | 400205

Protein expression 8, 18, 19,

Mutational profile P53

Culture Medium DMEM, w: 4.5g/L , w: 4mM L- , w: 3.7g/L NaHCO3, w: 1.0mM (820300a)

Supplements 10% FBS

Dissociation Reagent

Subculturing PBS . T25 3~5ml, T75 5~10ml PBS . T25

Fluid renewal 3~5

Freeze medium (FBS) + 10% DMSO , CM-1(Cytion 800100)

Thawing and Culturing Cells

1. .
2. -150°C , 3 .
3. 37°C 40~60 .
4. , 70% .
5. 8ml 15ml .
6. 300 x g 3 .
7. 10ml . T25 , T25
8. .

Hep-64.1 | 400205

Incubation Atmosphere 37°C, 5% CO₂

Flask Coating

Freezing Procedure -78°C

Shipping Conditions -78°C

Storage Conditions -150°C -196°C -80°C

/ / HLA

Sterility PCR