

NCI-H1792 Cells | 305835

General information

Description

NCI-H1792 is a human non-small cell lung carcinoma (NSCLC) cell line derived from a lung adenocarcinoma of an adult patient. It has been utilized extensively in cancer research, particularly in studies focusing on lung tumorigenesis, genetic aberrations, and drug sensitivity profiling. The cell line is characterized by an epithelial morphology and forms adherent monolayers in culture. Its inclusion in large-scale datasets such as the Cancer Cell Line Encyclopedia (CCLE) has enabled extensive genomic and proteomic profiling, facilitating comparative analyses with other lung cancer models.

Genomically, NCI-H1792 exhibits several molecular alterations common in NSCLC. It is known to harbor a KRAS mutation, a common oncogenic driver in lung adenocarcinoma, which contributes to aberrant MAPK signaling. The cell line has also been analyzed in proteomic studies, where its protein expression profile has provided insight into signaling pathway dependencies and vulnerabilities. Proteomic data highlight its utility in understanding pathway regulation and drug target validation across KRAS-mutant cancers. These datasets also underscore its classification within a subtype of KRAS-driven cancers that show distinct metabolic and signaling characteristics.

Organism

Human

Tissue

Metastatic

Disease

Lung adenocarcinoma

Synonyms

H1792, H-1792, NCIH1792

Characteristics

Age

50 years

Gender

Male

Ethnicity

Caucasian

Cell type

Epithelial

Growth properties

Adherent

Regulatory Data

Citation

NCI-H1792 (Cytion catalog number 305835)

## NCI-H1792 Cells | 305835

**Biosafety level** 1**NCBI\_TaxID** 9606**CellosaurusAccession** CVCL\_1495**Biomolecular Data****Mutational profile** Mutation: CDKN2A, Simple, p.Trp110Ter (c.330G>A) (p.Gly125Arg, c.373G>A), Heterozygous.Mutation, KRAS, Simple, p.Gly12Cys (c.34G>T), Heterozygous, TP53, Simple, c.672+1G>A, Homozygous, Note=Splice donor mutation**Handling****Culture Medium** RPMI 1640, w: 2.0 mM stable Glutamine, w: 2.0 g/L NaHCO<sub>3</sub> (Cytion article number 820700a)**Supplements** Supplement the medium with 10% FBS**Dissociation Reagent** Accutase**Doubling time** 45 hours**Fluid renewal** 2 to 3 times per week**Freeze medium** As a cryopreservation medium, we use complete growth medium (including FBS) + 10% DMSO for adequate post-thaw viability, or CM-1 (Cytion catalog number 800100), which includes optimized osmoprotectants and metabolic stabilizers to enhance recovery and reduce cryo-induced stress.

### NCI-H1792 Cells | 305835

#### Thawing and Culturing Cells

1. Confirm that the vial remains deeply frozen upon delivery, as cells are shipped on dry ice to maintain optimal temperatures during transit.
2. Upon receipt, either store the cryovial immediately at temperatures below -150°C to ensure the preservation of cellular integrity, or proceed to step 3 if immediate culturing is required.
3. For immediate culturing, swiftly thaw the vial by immersing it in a 37°C water bath with clean water and an antimicrobial agent, agitating gently for 40-60 seconds until a small ice clump remains.
4. Perform all subsequent steps under sterile conditions in a flow hood, disinfecting the cryovial with 70% ethanol before opening.
5. Carefully open the disinfected vial and transfer the cell suspension into a 15 ml centrifuge tube containing 8 ml of room-temperature culture medium, mixing gently.
6. Centrifuge the mixture at 300 x g for 3 minutes to separate the cells and carefully discard the supernatant containing residual freezing medium.
7. Gently resuspend the cell pellet in 10 ml of fresh culture medium. For adherent cells, divide the suspension between two T25 culture flasks; for suspension cultures, transfer all the medium into one T25 flask to promote effective cell interaction and growth.
8. Adhere to established subculture protocols for continued growth and maintenance of the cell line, ensuring reliable experimental outcomes.

#### Incubation Atmosphere

37°C, 5% CO<sub>2</sub>, humidified atmosphere.

#### Shipping Conditions

Cryopreserved cell lines are shipped on dry ice in validated, insulated packaging with sufficient refrigerant to maintain approximately -78 °C throughout transit. On receipt, inspect the container immediately and transfer vials without delay to appropriate storage.

#### Storage Conditions

For long-term preservation, place vials in vapor-phase liquid nitrogen at about -150 to -196 °C. Storage at -80 °C is acceptable only as a short interim step before transfer to liquid nitrogen.

## Quality Control & Molecular Analysis

**NCI-H1792 Cells | 305835**

**Sterility**

Mycoplasma contamination is excluded using both PCR-based assays and luminescence-based mycoplasma detection methods.

To ensure there is no bacterial, fungal, or yeast contamination, cell cultures are subjected to daily visual inspections.