

BC-3C Cells | 305246

General information

Description

BC3c is a human bladder carcinoma cell line established from a surgical biopsy of an invasive solid transitional cell carcinoma derived from an adult female donor. The cell line was generated from primary tumor tissue through enzymatic dissociation and serial passaging in McCoy's 5A medium supplemented with fetal calf serum. BC3c cells display an epithelial morphology with adherent growth and form cohesive colonies in vitro. Immunocytochemical characterization confirmed epithelial origin through expression of cytokeratins 8 and 19, while mesenchymal markers such as vimentin and fibroblast-associated antigens were absent. The cells do not express α -fetoprotein or carcinoembryonic antigen. Cytogenetic analysis demonstrated a near-triploid karyotype with multiple structural and numerical chromosomal abnormalities, consistent with genomic instability typical of invasive urothelial carcinoma.

Molecular analysis did not detect activating mutations in H-ras codon 12 or K-ras codons 12 and 13, and no mutations were identified in examined hotspot codons of TP53. Furthermore, there was no accumulation of p53 protein detected by immunocytochemistry, suggesting the absence of common p53 stabilization events in this model. BC3c cells exhibit rapid in vitro proliferation, including growth under low-serum conditions, indicative of autocrine growth stimulation. Conditioned medium from BC3c cultures supports proliferation, supporting the presence of endogenous growth-promoting factors.

Organism

Human

Tissue

Urinary bladder

Disease

Bladder carcinoma

Synonyms

BC3c

Characteristics

Age

82 years

Gender

Female

Ethnicity

Caucasian

Growth properties

Adherent

Regulatory Data

Citation

BC-3C (Cytion catalog number 305246)

BC-3C Cells | 305246

Biosafety level	1
NCBI_TaxID	9606
CellosaurusAccession	CVCL_1958

Biomolecular Data

Handling

Culture Medium	McCoy's 5a, w: 3.0 g/L Glucose, w: stable Glutamine, w: 2.0 mM Sodium pyruvate, w: 2.2 g/L NaHCO ₃ (Cytion article number 820200a)
Supplements	Supplement the medium with 10% FBS
Dissociation Reagent	Accutase 20 min 37°C
Doubling time	~25-30 hours
Freeze medium	As a cryopreservation medium, we use complete growth medium + 10% DMSO for adequate post-thaw viability.

BC-3C Cells | 305246

Thawing and Culturing Cells

1. Confirm that the vial remains deeply frozen upon delivery, as cells are shipped on dry ice to maintain optimal temperatures during transit.
2. Upon receipt, either store the cryovial immediately at temperatures below -150°C to ensure the preservation of cellular integrity, or proceed to step 3 if immediate culturing is required.
3. For immediate culturing, swiftly thaw the vial by immersing it in a 37°C water bath with clean water and an antimicrobial agent, agitating gently for 40-60 seconds until a small ice clump remains.
4. Perform all subsequent steps under sterile conditions in a flow hood, disinfecting the cryovial with 70% ethanol before opening.
5. Carefully open the disinfected vial and transfer the cell suspension into a 15 ml centrifuge tube containing 8 ml of room-temperature culture medium, mixing gently.
6. Centrifuge the mixture at 200 x g for 5 minutes, carefully discard the supernatant containing freezing medium.
7. Follow the procedure described under Post-Thaw Recovery

Incubation Atmosphere

37°C, 5% CO₂, humidified atmosphere.

Shipping Conditions

Cryopreserved cell lines are shipped on dry ice in validated, insulated packaging with sufficient refrigerant to maintain approximately -78 °C throughout transit. On receipt, inspect the container immediately and transfer vials without delay to appropriate storage.

Storage Conditions

For long-term preservation, place vials in vapor-phase liquid nitrogen at about -150 to -196 °C. Storage at -80 °C is acceptable only as a short interim step before transfer to liquid nitrogen.

Quality Control & Molecular Analysis