

FBS, Fetal Bovine Serum, SA-Origin, 500 ml | 831109

Sourced from South American bovine sources, this sterile-filtered Fetal Bovine Serum (FBS) provides consistent cell culture performance. Prior to release, every batch is tested for sterility, mycoplasma, and bovine viruses. Furthermore, each lot is evaluated to verify its ability to support the growth of HEK293 and CHO cells. The serum is available in 50 mL and 500 mL volumes, and a lot-specific Certificate of Analysis (CoA) can be provided upon request.

Product Specifications

Appearance	Clear amber liquid
Storage and shelf life	Store at $\leq -15^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles. Preparation of aliquots recommended. Once opened store at 4°C and use within 4-6 weeks.
Shipping conditions	Frozen (Dry ice)
Thawing	Overnight at $+2^{\circ}\text{C}$ to $+8^{\circ}\text{C}$. Swirl gently to homogenize.