

FBS, Fetal Bovine Serum, US-Origin, 50 ml | 831108

Sourced from the United States, this sterile-filtered Fetal Bovine Serum (FBS) provides a consistent foundation for cell culture applications. Every batch undergoes quality control screening for sterility, mycoplasma, and bovine viruses. Additionally, each lot is biologically validated to confirm its ability to support the growth of CHO and HEK293 cell lines. The serum is packaged in 50 mL and 500 mL volumes. Batch-specific Certificates of Analysis (CoA) are available upon request.

Product Specifications	
Appearance	Clear amber liquid
Storage and shelf life	Store at $\leq -15^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles. Preparation of aliquots recommended. Once opened store at 4°C and use within 4-6 weeks.
Shipping conditions	Frozen (Dry ice)
Thawing	Overnight at $+2^{\circ}\text{C}$ to $+8^{\circ}\text{C}$. Swirl gently to homogenize.