

MC3T3-E1 Subclone 24 Cells | 305186**General information****Description**

MC3T3-E1 Subclone 24 cells expressly represent a preosteoblast cell type, which plays a crucial role in bone formation. Morphologically, they exhibit a fibroblast-like appearance, characterized by their elongated shape and spindle-like structures. This particular subclone is derived from the calvaria tissue, a skull region contributing to bone formation. One of the critical applications of MC3T3-E1 Subclone 24 Cells lies in 3D cell culture, where researchers can study the behavior and interactions of these cells within a three-dimensional environment. This method offers a more physiologically relevant model than traditional two-dimensional cell cultures, allowing for a better understanding of the intricate processes involved in bone formation.

While these cells possess numerous advantages, it's important to note their specific characteristics. MC3T3-E1 Subclone 24 Cells have been observed to exhibit poor osteoblast differentiation when exposed to ascorbic acid, a key component for promoting bone cell growth. Furthermore, they do not form a mineralized extracellular matrix, a crucial step in creating bone tissue. The doubling time of MC3T3-E1 Subclone 24 Cells is approximately 90.5 hours.

Organism Mouse**Tissue** Bone**Applications** 3D cell culture, Differentiation studies**Characteristics****Breed/Subspecies** C57BL/6**Age** 1 day**Gender** Unspecified**Morphology** Fibroblast**Cell type** Osteoblast**Growth properties** Adherent**Regulatory Data****Citation** MC3T3-E1 Subclone 24 (Cytion catalog number 305186)**Biosafety level** 1

MC3T3-E1 Subclone 24 Cells | 305186**NCBI_TaxID** 10090**CellosaurusAccession** CVCL_5438**Biomolecular Data****Receptors expressed** Parathyroid hormone-related protein (PTHrP) receptor**Protein expression** Collagen, bone sialoprotein (BSP), osteocalcin (OCN), parathyroid hormone (PTH)**Tumorigenic** Yes, in immunosuppressed mice**Handling****Culture Medium** Alpha MEM, w: 2.0 mM stable Glutamine, w: Ribonucleosides, w: Deoxyribonucleosides, w: 1.0 mM Sodium pyruvate, w: 2.2g/L NaHCO₃, w/o: Ascorbic acid (GIBCO, Catalog No. A1049001. We do not supply this product; please consider other suppliers. Please let us know if you need further assistance.)**Supplements** Supplement the medium with 10% FBS**Dissociation Reagent** Accutase**Subculturing** Remove the old medium from the adherent cells and wash them with PBS that lacks calcium and magnesium. For T25 flasks, use 3-5 ml of PBS, and for T75 flasks, use 5-10 ml. Then, cover the cells completely with Accutase, using 1-2 ml for T25 flasks and 2.5 ml for T75 flasks. Let the cells incubate at room temperature for 8-10 minutes to detach them. After incubation, gently mix the cells with 10 ml of medium to resuspend them, then centrifuge at 300xg for 3 minutes. Discard the supernatant, resuspend the cells in fresh medium, and transfer them into new flasks that already contain fresh medium.**Freeze medium** As a cryopreservation medium, we use complete growth medium (including FBS) + 10% DMSO for adequate post-thaw viability, or CM-1 (Cytion catalog number 800100), which includes optimized osmoprotectants and metabolic stabilizers to enhance recovery and reduce cryo-induced stress.

MC3T3-E1 Subclone 24 Cells | 305186

Thawing and Culturing Cells

1. Confirm that the vial remains deeply frozen upon delivery, as cells are shipped on dry ice to maintain optimal temperatures during transit.
2. Upon receipt, either store the cryovial immediately at temperatures below -150°C to ensure the preservation of cellular integrity, or proceed to step 3 if immediate culturing is required.
3. For immediate culturing, swiftly thaw the vial by immersing it in a 37°C water bath with clean water and an antimicrobial agent, agitating gently for 40-60 seconds until a small ice clump remains.
4. Perform all subsequent steps under sterile conditions in a flow hood, disinfecting the cryovial with 70% ethanol before opening.
5. Carefully open the disinfected vial and transfer the cell suspension into a 15 ml centrifuge tube containing 8 ml of room-temperature culture medium, mixing gently.
6. Centrifuge the mixture at $300 \times g$ for 3 minutes to separate the cells and carefully discard the supernatant containing residual freezing medium.
7. Gently resuspend the cell pellet in 10 ml of fresh culture medium. For adherent cells, divide the suspension between two T25 culture flasks; for suspension cultures, transfer all the medium into one T25 flask to promote effective cell interaction and growth.
8. Adhere to established subculture protocols for continued growth and maintenance of the cell line, ensuring reliable experimental outcomes.

Incubation Atmosphere

37°C , 5% CO_2 , humidified atmosphere.

Flask Coating

None

Shipping Conditions

Cryopreserved cell lines are shipped on dry ice in validated, insulated packaging with sufficient refrigerant to maintain approximately -78°C throughout transit. On receipt, inspect the container immediately and transfer vials without delay to appropriate storage.

Storage Conditions

For long-term preservation, place vials in vapor-phase liquid nitrogen at about -150 to -196°C . Storage at -80°C is acceptable only as a short interim step before transfer to liquid nitrogen.

Quality Control & Molecular Analysis

MC3T3-E1 Subclone 24 Cells | 305186

Sterility

Mycoplasma contamination is excluded using both PCR-based assays and luminescence-based mycoplasma detection methods.

To ensure there is no bacterial, fungal, or yeast contamination, cell cultures are subjected to daily visual inspections.