

NCI-H2009 | 305283

Thawing and Culturing Cells

1. [REDACTED]
2. [REDACTED] cryovial [REDACTED]-150 [REDACTED]
3. [REDACTED] 37°C [REDACTED]
4. [REDACTED] 70%
5. [REDACTED] 15 [REDACTED]
6. [REDACTED] 300 x g [REDACTED] 3 [REDACTED]
7. [REDACTED] 10 [REDACTED]
8. [REDACTED]

Incubation Atmosphere 37°C, 5% CO₂, [REDACTED]

Flask Coating None

Shipping Conditions [REDACTED]

Storage Conditions [REDACTED]-150 [REDACTED]-196 °C [REDACTED]

[REDACTED] / [REDACTED] / HLA

Sterility [REDACTED] PCR [REDACTED]
 [REDACTED] [REDACTED] [REDACTED]