

Product sheet

MOLM-13 | 305393

Antigen expression CD3-, CD4+, CD14-, CD15+, CD19-, CD33+, CD34-, cy CD68+, HLA-DR-

Mutational profile FLT3, tandem; KMT2A-MLLT3, MLL-MLLT3, MLL-AF9

Culture Medium RPMI 1640, w: 2.0 mM, w: 2.0 NaHCO₃ (Cytion 820700a)

Supplements FBS 10%

Seeding density 5 6

Fluid renewal 2 3

Freeze medium (FBS) + 10% DMSO

- Thawing and Culturing Cells**
- 1.
 2. cryovial -150
 3. 37°C
 4. 70%
 5. 15
 6. 300 x g 3
 7. 10
 - 8.

Incubation Atmosphere 37°C, 5% CO₂

