

|                      |                                |
|----------------------|--------------------------------|
| Description          | MV3 (xenograft models) MV3 MV3 |
| Organism             |                                |
| Tissue               |                                |
| Disease              |                                |
| Metastatic site      |                                |
| Synonyms             | MV-3, MV 3                     |
| Age                  | 76                             |
| Gender               |                                |
| Ethnicity            |                                |
| Morphology           |                                |
| Growth properties    |                                |
| Citation             | MV3 ( Cytion 305495)           |
| Biosafety level      | 1                              |
| NCBI_TaxID           | 9606                           |
| CellosaurusAccession | CVCL_W280                      |

MV3 | 305495

|                             |                                                                                                                     |
|-----------------------------|---------------------------------------------------------------------------------------------------------------------|
| Culture Medium              | DMEM, w: 4.5 g/l / 4.5 g/l, w: 4 g/l - 4 g/l, w: 3.7 g/l / 3.7 g/l NaHCO <sub>3</sub> , w: 1.0 g/l                  |
| Supplements                 | FBS 10%                                                                                                             |
| Dissociation Reagent        |                                                                                                                     |
| Subculturing                | PBS T25                                                                                                             |
| Freeze medium               | (FBS) + 10% DMSO                                                                                                    |
| Thawing and Culturing Cells | <div><div>1.</div><div>2.</div><div>3.</div><div>4.</div><div>5.</div><div>6.</div><div>7.</div><div>8.</div></div> |
| Incubation Atmosphere       | 37°C, 5% CO <sub>2</sub>                                                                                            |
| Shipping Conditions         |                                                                                                                     |

MV3 | 305495

Storage  
Conditions

Store at -150 °C to -196 °C

Sterility

PCR ready  
No DNase treatment required