

XXXXXXXXXXXXXXXXXXXX

XXXXXXXXXX

XX

CellosaurusAccession CVCL_0626

[illegible]XXXXXXXXXXXX

Cytion GmbH | XXXXXXXX-XXXXX3 | 69124 XXXXXXXXXX | XXXXXXXX
 XXXX +49(0)6221 405780 | www.cytion.com | info@cytion.com

SK-CO-1 | 305626

Thawing and
Culturing Cells

- 1. Thaw cells rapidly at 37°C.
- 2. Add cells to 5 ml cryovial containing 1 ml of medium.
- 3. Incubate cells at 37°C for 24 hours.
- 4. Seed cells into 24-well plate at 70% confluency.
- 5. Incubate cells at 37°C for 24 hours.
- 6. Harvest cells by centrifugation at 300 x g for 3 minutes.
- 7. Wash cells with PBS.
- 8. Resuspend cells in 100 µl of medium.

Incubation
Atmosphere

37°C, 5% CO₂, humidified

Shipping
Conditions

Cells are shipped on dry ice in a cool box.

Storage
Conditions

Cells can be stored at -150 °C or -196 °C in liquid nitrogen.

Sterility

Cells are tested for mycoplasma contamination by PCR. Cells are also tested for endotoxin levels.