

XXXXXXXXXXXXXXXXXXXX

Z-138 (MCL) BCL-2 PI3K/AKT

XXXXXXXXXX

Z138

70 ㊸㊸



XXXXXXXXXX

XXXXXXXXXXXXXXXXXXXX

[illegible]

XXXXXXXXXXXXXXXXXXXXXXXXXXXX

Z-138 (XXXXXXXXXXXXXXXXXXXX Cytion 305516)

1

9606

Product sheet

XXXXXXXXZ-138 | 305516

CellosaurusAccession CVCL_B077

XXXXXXXXXXXXXXXXXXXX

Antigen expression CD3-, CD5-, CD10-, CD19+, CD20+, CD23+, FMC7-

Mutational profile XXXXXXXXXXXX TRAF2, XXXXXXXXXXXX p.Trp114Ter (c.341G>A), XXXXXXXX

XXXXXXXXXXXX

Culture Medium IMDM, w: 4.5 XXXX/XXXX XXXXXXX, w: 4 XXXXXXXXXXXX XXX-XXXXXXX, w: 25 XXXXXXXXXXXX HEPES, w: 1.0 XXXXXXXXXXXX XXXXXXXXXXXX

Supplements XXXXXXXXXXXX 10% XXXXXXXXXXXXXXXXXXXXXXX

Doubling time 24 XXXXXXX

Subculturing XXXXX⁵ XXXXX⁵ XXXXX⁶ XXX 5 x 10

Seeding density 3–5 × 10⁵/XXXXXXXX

Fluid renewal 2 XXX 3 XXXXXXXXXXXXXXX

Freeze medium XXX (XXXXXX FBS) + 10% DMSO

XXXXXX Z-138 | 305516


Thawing and Culturing Cells

1.              

Incubation Atmosphere

37°C, 5% CO₂,

Shipping Conditions



Storage Conditions

-150 -196 °C

Sterility

Age Group	PCR (%)
0-9	5
10-19	10
20-29	85
30-39	75
40-49	65
50-59	55
60-69	45
70-79	35
80-89	25
90-99	15

[illegible]